




KolXPD

software for photoelectron spectroscopy
data measurement and processing

FEATURE HIGHLIGHTS

- Focused to photoelectron spectroscopy
- Very fast data viewing and processing
- Energy- and angle-resolved experiments
- Import and export of variety formats
- Huge amount of data handling
- Fitting, waterfall charts, quantitative preview
- Functionality customizable on request



Charles University
in Prague, Czech Republic



Academy of Sciences of the
Czech Republic, Prague



Synchrotron Elettra
Trieste, Italy



National Institute
for Materials Science
Tsukuba, Japan



Synchrotron SPring-8
Sayo, Japan




Heinrich-Heine-Universität
Düsseldorf, Germany

most significant institutes where KolXPD is used


Photoelectron spectroscopy features

- import of data from many analyzer softwares (SpecsLab, Spectra, Omicron, VG, etc.)
- complete XML experiment file of SpecsLab can be loaded and processed, including each region channel and sweep data
- spectra can be acquired by KolXPD (currently with Specs analyzer, others on request)
- support for synchrotron radiation spectroscopy
- spectra can be quickly fitted with most common peak shapes and backgrounds
- all data are organized in one tree structure
- just select more items to view them in one multiview chart
- multiview has many icons for instant processing which doesn't affect original data. Just click to view normalized spectra, make waterfall, subtract background, show fit results, smooth, etc.
- using multiview it is possible to fit dozens of spectra quickly and automatically
- all result data can be exported, all graphs copied into clipboard and pasted anywhere
- data can be transferred in both directions between KolXPD and WaveMetrics Igor Pro for further processing and presentation





Multidimensional experiments

- support for energy-resolved experiments using synchrotron radiation – methods like NEXAFS, Resonant photoemission, etc.
- support for angle-resolved spectroscopy, e.g. XPD (X-ray Photoelectron Diffraction) and ARUPS (Angle-Resolved UV Photoelectron Spectroscopy - band mapping)
- KolXPD can measure these experiments using analyzer, sample manipulator or beamline monochromator
- 2-dimensional angle-resolved data can be plotted as a pattern picture in Cartesian or polar coordinates
- thousands of 2D experiment spectra can be easily and automatically fitted by one shape
- intensity of each pattern point can be spectrum area or any parameter of fit function
- results of XPD pattern modeling software EDAC can be imported and displayed in same way
- multidimensional data can be imported and exported, including pattern pictures – as both text and image format



Other features

- tree structure is useful for whole experiment handle; other items can be added to describe experiment steps e.g. cleaning, heating, deposition, gas exposure etc.
- you can add own notes to any item; folders and note items can contain HTML-formatted text
- simple image browser e.g. for LEED pictures
- KolXPD allows simple quantitative preview by showing peak areas corrected by analyzer transmission function and photoemission cross-section
- relative concentrations can be showed for each folder; it's useful for preview of concentration evolution during experimental steps
- all files of KolXPD has open text format – you can write an import to your favorite software



Watch video-tutorial at

<http://www.kolibrik.net/kolxpd/>


How to import 15 spectra, make graphs, fit peaks and analyze results in 90 seconds !!!

Tutorial steps:


- start: 15 raw XY text files with sample spectra (2 doublets)
- import to KolXPD by single drag-and-drop
- place spectra to one graph, change x-axis to binding energy
- change display to waterfall and color-filled peak area
- define x-axis range and fitting function for middle spectrum
- automatically fit all spectra with this function
- display fit results
- export area evolution of doublets and display in new graph
- fit area evolution with linear function

All this done in 90 seconds!


Video snapshots:




Raw XY text files



Waterfall view in KolXPD






Fitting of one spectrum



Doublets - area evolution


Angle-resolved photoemission experiments acquired using KolXPD at Charles University in Prague

2D angle-resolved spectra of core-levels and valence band can be measured. All measured spectra can be automatically fitted and components of fit can be used for displaying a pattern. Valence band can be displayed as polar or azimuthal plot or as 2D map like Fermi surface map.





XPD pattern calculations made by EDAC* patterns and clusters can be visualized by KolXPD

*Software for X-ray photoelectron diffraction (XPD) pattern calculations
F. J. García de Abajo, M. A. Van Hove, and C. S. Fadley, Phys. Rev. B 63, 075404 (2001).




Angle- and energy-resolved photoemission experiments acquired using KolXPD at Materials Science Beamline, synchrotron Elettra

KolXPD can control both manipulator and monochromator – energy-resolved experiments like NEXAFS, resonant photoemission, etc. are possible. Graphs below show resonant experiment on cerium-oxide and XPD patterns acquired at two different resonant photon energies.



Reference papers and results where KolXPD was used for measurement and data processing

- **J. Libra, V. Matolín, Surf. Sci. 600 (2006) 2317**
Angle resolved photoemission study of the Ce/Pd(111) interface
- **J. Libra, et al., Phys. Rev. B 76 (2007) 165438**
Angle-resolved photoemission study of ordered Pb/Ni(111) surface phases
- **V. Matolín, J. Libra, Surf. Sci. 601 (2007) 4058–4062**
Fermi surface and band mapping of the cerium/palladium surface alloy
- **M. Škoda, J. Libra, et al., Surface Science 601 (2007) 4958**
A resonant photoemission study of the Ce and Ce-oxide/Pd(111) interfaces
- **V. Matolín, J. Libra, et al., J. Phys. Chem. C 112 (2008) 3751**
Photoemission spectroscopy study of Cu/CeO₂ systems: Cu/CeO₂ nanosized catalyst and CeO₂(111)/Cu(111) inverse model catalyst



(sample images made using results of KolXPD data processing)

About author

Mgr. Jiří Libra, Ph.D. (email: Jiri.Libra@gmail.com)

- Surface science experienced during 10 years spent in Department of Surface and Plasma Science at Faculty of Mathematics and Physics, Charles University in Prague under supervision of *prof. Vladimír Matolín*
- Business activities: development of software and electronics (not only) for science. Experiment automation, data acquisition and processing, electronics for angle-resolved experiments, software for beamline monochromator control, scientific device interfacing, etc.

KolXPD Web

<http://www.kolibrik.net/kolxpd/>

- You can download KolXPD and track new versions, download documents, manuals, tutorials and this poster
- Purchase a license or make a donation, get a trial license for evaluation
- Contact author, visit KolXPD forum, discuss with other KolXPD users